

LINGUAL INDIRECT BONDING / Rx



<input type="checkbox"/> PHONE ME REGARDING THIS CASE	SEND ADDITIONAL <input type="checkbox"/> RX SHEETS <input type="checkbox"/> MAILING LABELS
<input type="checkbox"/> SPECIAL INSTRUCTIONS ON FILE	
<input type="checkbox"/> NEW ACCOUNT <input type="checkbox"/> ADDRESS CHANGE	

Doctor _____

Address _____

City _____ State _____ Zip _____


Telephone _____ Fax _____

Patient Name _____

Date Shipped _____ Date Needed _____

Appointment Time _____

e-Mail Address _____



CASE INFORMATION	Upper	Lower	Slot Size
LINGUAL	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> .018
TRAY SECTIONS	<input type="checkbox"/> 1 Piece	<input type="checkbox"/> Midline Split	<input type="checkbox"/> 3 Pieces / Arch

Please Indicate on Diagram to Right

1. Circle teeth to be banded
2. Mark an "X" on teeth missing, to be extracted, or those not to be bonded
3. Indicate with arrows over-rotations desired

BRACKET PLACEMENT INFORMATION - ALL CASES

Please use chart to indicate changes from standard values

	STANDARD VALUES			CHANGES		
	TQ	ANG	HT	TQ	ANG	HT
<u>1 1</u>	14	5	X	_____	_____	_____
<u>2 2</u>	7	8	X-.5mm	_____	_____	_____
<u>3 3</u>	-2	10	X+.5mm	_____	_____	_____
<u>54 45</u>	-7	0	X	_____	_____	_____
<u>76 67</u>	-10	0	X	_____	_____	_____
<u>21 12</u>	0	0	X	_____	_____	_____
<u>3 3</u>	-7	6	X+.5mm	_____	_____	_____
<u>4 4</u>	-12	0	X	_____	_____	_____
<u>5 5</u>	-16	0	X	_____	_____	_____
<u>6 6</u>	-22	0	X	_____	_____	_____
<u>7 7</u>	-27	0	X	_____	_____	_____

LINGUAL CASE INFORMATION

Lingual Bracket Placement Method

3x3 CLASS Set-Up TARG

Lingual Archwires Requested

Indicate on Reverse Side

Ormco Lingual Molar Codes - Please Mark Appropriate Molar Brackets on Above Diagram

TT Terminal Tube **TB** Twin Bracket
HC Hinge Cap **BTS** Twin w/ TP Sheath

Other Lingual Services

Specialty Supply and Fit Molar Bands
 Specialty Weld Molar Attachments to Bands
 Soldered Transpalatal Wire to Molars
 Removable Transpalatal Wire to Molars

SPECIAL INSTRUCTIONS _____



800-423-3270
www.ApplianceTherapy.com